

ATOMIZER



GUARANTEED
3
YEARS

AEROSOL CON COMPRESSORE A PISTONE PER USO ORDINARIO PISTON COMPRESSOR NEBULIZER FOR HOME USE

DATI TECNICI COMPRESSORE 2 VELOCITÀ

COMPRESSORE	A PISTONE
PRESSIONE MASSIMA	185 kPa
PRESSIONE D'ESERCIZIO	I VEL. 45 kPa - II VEL. 60 kPa
PORTATA D'ARIA	I VEL. 8 l/min - II VEL. 11 l/min
CICLO DI UTILIZZO	CONTINUO
RUMOROSITÀ	56 dBA
PESO	1,4 kg
CLASSE DI RISCHIO	IIa

ACCESSORI

AMPOLLA NEBJET (POLIPROPILENE)
CONNETTORE PER ACCESSORI (POLIPROPILENE)
NASALE (POLIPROPILENE)
BOCCAGLIO (POLIPROPILENE)
MASCHERA ADULTO (PVC DOP FREE)
MASCHERA PEDIATRICA (PVC DOP FREE)
TUBO ARIA (PVC DOP FREE)
FILTRI DI RICAMBIO (3 PZ.)

COMPRESSOR TECHNICAL DATA 2 SPEEDS

COMPRESSOR	PISTON COMPRESSOR
MAX PRESSURE	185 kPa
OPERATING PRESSURE	I SPEED 45 kPa - II SPEED 60 kPa
AIR FLOW	I SPEED 8 l/min - II SPEED 11 l/min
OPERATING CYCLE	CONTINUOUS
NOISE LEVEL	56 dBA
WEIGHT	1.4 kg
MDD CLASS	IIa

ACCESSORIES

NEBJET NEBULIZER (POLYPROPYLENE)
ACCESSORIES'S CONNECTOR (POLYPROPYLENE)
NOSEPIECE (POLYPROPYLENE)
MOUTHPIECE (POLYPROPYLENE)
ADULT MASK (PVC DEPH FREE)
CHILD MASK (PVC DEPH FREE)
AIR TUBE (PVC DEPH FREE)
SPARE FILTERS (3 PCS.)

C€0434

(dir. 93/42/EEC)

PCS.

6

MASTER CARTON



DATI TECNICI AMPOLLA

CAPACITÀ	12 ml
FRAZIONE RESPIRABILE	< 5 μ m 75%
MMAD (Mass Median Aerodynamic Diameter)	2,61 μ m
NEB RATE I vel.	0,20 ml/min
NEB RATE II vel.	0,35 ml/min
VOLUME RESIDUO	0,7 ml

NEBULIZER TECHNICAL DATA

CAPACITY	12 ml
RESPIRABLE VOLUME	< 5 μ m 75%
MMAD (Mass Median Aerodynamic Diameter)	2.61 μ m
NEB RATE I SPEED	0.20 ml/min
NEB RATE II SPEED	0.35 ml/min
RESIDUAL VOLUME	0.7 ml

