



GUARANTEED
3
YEARS

**AEROSOL CON COMPRESSORE A PISTONE PER USO ORDINARIO
PISTON COMPRESSOR NEBULIZER FOR HOME USE**

DATI TECNICI COMPRESSORE

COMPRESSORE	A PISTONE
PRESSIONE MASSIMA	185 kPa
PRESSIONE D'ESERCIZIO	60 kPa
PORTATA D'ARIA	10 l/min
CICLO DI UTILIZZO	20/40 min. I/O
RUMOROSITÀ	58 dBA
PESO	1,1 kg
CLASSE DI RISCHIO	Ila

COMPRESSOR TECHNICAL DATA

COMPRESSOR	PISTON COMPRESSOR
MAX PRESSURE	185 kPa
OPERATING PRESSURE	60 kPa
AIR FLOW	10 l/min
OPERATING CYCLE	20/40 min. I/O
NOISE LEVEL	58 dBA
WEIGHT	1.1 kg
MDD CLASS	Ila



MASCHERA ADULTO - ADULT MASK



TUBO ARIA - AIR TUBE



AMPOLLA QUICKJET - QUICKJET NEBULIZER

CE0434
(dir. 93/42/EEC)

PCS. **6**
MASTER CARTON



MASCHERA PEDIATRICA - CHILD MASK



FILTRI - FILTERS

BOCCAGLIO / NASALE INTEGRATI - MOUTHPIECE / NOSEPIECE INTEGRATED



2 IN 1



RACCORDO ACCESSORI
ACCESSORIE'S CONNECTOR

ACCESSORI

- AMPOLLA QUICKJET (POLIPROPILENE)
- CONNETTORE PER ACCESSORI (POLIPROPILENE)
- BOCCAGLIO / NASALE INTEGRATI (POLIPROPILENE)
- MASCHERA ADULTO (PVC DOP FREE)
- MASCHERA PEDIATRICA (PVC DOP FREE)
- TUBO ARIA (PVC DOP FREE)
- FILTRI DI RICAMBIO (3 PZ.)

ACCESSORIES

- QUICKJET NEBULIZER (POLYPROPYLENE)
- ACCESSORIE'S CONNECTOR (POLYPROPYLENE)
- NOSEPIECE / MOUTHPIECE INTEGRATED (POLYPROPYLENE)
- ADULT MASK (PVC DEPH FREE)
- CHILD MASK (PVC DEPH FREE)
- AIR TUBE (PVC DEPH FREE)
- SPARE FILTERS (3 PCS.)

DATI TECNICI AMPOLLA

CAPACITÀ	9 ml
FRAZIONE RESPIRABILE	< 5 µm 75%
MMAD (Mass Median Aerodynamic Diameter)	2,61 µm
NEB RATE	0,35 ml/min
VOLUME RESIDUO	0,7 ml

NEBULIZER TECHNICAL DATA

CAPACITY	9 ml
RESPIRABLE VOLUME	< 5 µm 75%
MMAD (Mass Median Aerodynamic Diameter)	2.61 µm
NEB RATE	0.35 ml/min
RESIDUAL VOLUME	0.7 ml